

KEY ACCESSORIES

REFRIGERATED CHAMBERS

Models: 7001-10/25/33/50/75



One-Handed Door Latch

Secure yet easy chamber access. Pull handle to open, door latches automatically. Available with or without integral lock, for 25 & 33 ft³ units.



Extended Temperature Range

Extends unit operation down to 0°C, includes defrost cycle for 0-4.5°C



Additional Shelving

Perforated, solid, wire, and heavy-duty variants. All shelves made of polished type 304 stainless steel and include shelf brackets. Factory-installed kits available to swap all shelves for another type.



Outputs / Recorders

Circular temp/RH chart recorders in 6", 10" and 12" sizes, plus analog outputs and digital USB datalogging.