

Q1 2008 NEWSLETTER

MARCH 31, 2008

VOLUME 6, NUMBER 1

CARON APPOINTS EUROPEAN SALES CONSULTANT

Caron Products has partnered with Bob Batty International based in Cambridge, UK, who specialize in assisting distributors of laboratory products on behalf of select laboratory equipment manufacturers. BBI is focused on assisting current European Caron distributors, as well as finding new distributors in countries where we are not represented.

Mr. Batty will be representing Caron at Analytica 2008 in Munich, Germany in early April. He will be conducting meetings with existing distributors and will be interviewing potential companies.

Mr. Batty can be contacted at www.bobbatty.com or by email at international@bobbatty.com.

MULTIPLE OPTIONS AVAILABLE TO SATISFY SPECIFIC CUSTOMER NEEDS

In an effort to better serve our chamber customers, Caron has introduced the following accessory items:

- Cat. No. SHKR101 Stainless Steel Triple-tier Shaker Support System. The system is capable of supporting up to three 150 lb. shakers operating at 150 RPM each. There is 19" clearance between each shelf, and the Shaker Support System comes complete with two 115V GFI duplex outlets. When ordering a 30 cu. ft. chamber with this accessory, the standard leveling feet on the chambers are replaced with vibration isolation pads.



- Cat. No. CART102 Heavy-duty Laboratory Equipment Cart. This cart is 48" x 30" x 24" and has a reinforced, solid 12 gauge, removable aluminum top. The under shelf provides convenient storage space. The cart has 5" stem-type casters with brakes and can accommodate up to 800 lbs.—ideal for Caron's 10 cu. ft. sized chambers.

- Cat. No. TABL101 Stainless Steel Work Table. The fully welded 30" x 48" x 35" cabinet base has sliding doors, a 1-1/4" turn down lip on four sides, 14 gauge stainless steel top, adjustable legs, and under shelf. This unit ships fully assembled and is ideal for Caron's 10 cu. ft. sized chambers.



CARON Products & Services, Inc.
27640 St. Rt. 7
PO Box 715
Marietta, OH 45750

Phone
740-373-6809
800-648-3042

Fax
740-374-3760

E-mail
sales@caronproducts.com

We're on the Web!
www.caronproducts.com

Contact Caron for pricing & availability at (1) 740-373-6809 or by email at: sales@caronproducts.com

CARON
Opening Doors for Scientists

If you would prefer not to receive further messages from this sender:

• Click on the Reply button. • Replace the Subject field with the word REMOVE. • Click the Send button.
Thank You.