

NEWSLETTER

FEBRUARY 2009

VOLUME 7, NUMBER 1

A Family of Exciting New Products



CARON's all NEW family of chambers delivers more than accuracy, dependability and high performance. They're loaded with earth friendly technical innovations that provide precise temperature and humidity controlled conditions for numerous applications. The 25 cu. ft. units come complete with factory-installed, heavy-duty casters and will roll through a standard 34" x 80" doorway.

Environmental Chambers

These chambers were designed in accordance with ICH Guidelines (Q1A) and are ideal for stability and package testing, temperature and humidity elevation studies, shelf life studies, and more.

- Our new **gVapor™** humidity system delivers controlled humidity vapor on an as-needed basis without wasting energy or generating heat.
- **gROD™** (Refrigeration on Demand) control system efficiently manages power consumption and saves electricity.
- Units offer a wide temperature range of 5°C to 70°C and humidity range from 20% to 98%.
- Unique "tool-less" removal of interior components minimizes downtime for cleaning.
- The ergonomically-designed door handle allows easy access into the chamber.
- The 10 cu. ft. units offer unique "stackability" feature.



Model 6020



(2) Model 6010s

Reach-in CO₂ Incubator

Designed with customer convenience in mind, this unit features 25 cu. ft. of precise temperature, IR-CO₂, and controlled humidity space.

- Industry's first large capacity incubator to offer an automatic moist heat decontamination cycle.
- An IR (infrared) sensor is utilized to maintain a CO₂ range of 0 to 20%.
- The energy efficient **gVapor™** humidity system eliminates the need for a water pan and reduces contamination.
- With "tool-less" removal of interior components, users have virtually unlimited shelving configurations that are easily field installed to accommodate shakers, cell rollers, and other requirements.



Model 6024

CARON Welcomes New Director of Service & Quality



Jim Lang has joined the company as Director of Service & Quality. Lang brings 31 years of extensive laboratory equipment experience in both service and quality to CARON. In his new position, Jim is responsible for establishing a nationwide service network to assist CARON with its rapidly growing customer base. Additionally, he will lead CARON's efforts to implement more extensive quality control techniques in manufacturing.

CARON Products & Services, Inc.
27640 St. Rt. 7
PO Box 715
Marietta, OH 45750

Phone
740-373-6809
800-648-3042

Fax
740-374-3760

E-mail
sales@caronproducts.com

We're on the Web!
www.caronproducts.com

CARON
Opening Doors for Scientists

If you would prefer not to receive further messages from this sender:
• Click on the Reply button. • Replace the Subject field with the word REMOVE. • Click the Send button.
Thank You.