

KEY ACCESSORIES

INSECT REARING CHAMBERS

Models: 7340-25/33



One-Handed Door Latch

Secure yet easy chamber access. Pull handle to open, door latches automatically. Available with or without integral lock, for 25 & 33 ft³ units.



Hydrogen Peroxide Sterilization

Unique validated 2-hour dry H₂O₂ system. Provides 12 log sterilization with no post-cycle cleanup. Consists of H₂O₂ module, system prep, and disposable H₂O₂ canisters.



Controlled Humidity

Precise RH control without messy water pans or flooded unit bottoms. Simple atomizer-based & sensor-controlled system generates no heat and is low maintenance.



"No-Blue" Colored Lighting

Narrow 500-700 nanometer lighting eliminates lethal blue light effects documented in flying insects. LED modules are wall-mounted, sealed, and removable without tools for service or cleaning.



Full Spectrum Lighting

Low intensity white diurnal lighting. Efficient LED technology for long life and low power use. Wall-mounted sealed modules removable without tools for service or cleaning.



Condensate Recirculator

Provides a source of Type II equipment water and eliminates the need for a drain line. Includes sediment & DI filters, and a 7-gallon reservoir. Compact design fits between chambers and rolls out for easy access.