

# KEY ACCESSORIES

## HEATED/HUMIDIFIED INCUBATORS

*Models: 7002-10/25/33, 7003-10/25/33*



### One-Handed Door Latch

---

Secure yet easy chamber access. Pull handle to open, door latches automatically. Available with or without integral lock, for 25 & 33 ft<sup>3</sup> units.



### Column Air Dryer

---

Permits chamber operation down to very low RH setpoints. Requires external compressed dry air source. Noise suppression kit also available for in-lab use.



### Condensate Recirculator

---

Provides a source of Type II equipment water and eliminates the need for a drain line. Includes sediment & DI filters, and a 7-gallon reservoir. Compact design fits between chambers and rolls out for easy access.



### Additional Shelving

---

Perforated, solid, wire, and heavy-duty variants. All shelves made of polished type 304 stainless steel and include shelf brackets. Factory-installed kits available to swap all shelves for another type.