

BATH CIRCULATOR

Model: 2050

Model	2050-1	2050-4
Temperature Range	-15°C to 90°C	
Temperature Control	± 0.1°C	
Temperature Controller	Digital display, PID control, 1/16 DIN	
Cooling Capacity @ 20°C	1025 BTU/Hr, 300 Watts	
External Pumping	6.5 LPM at 0' head	
Work Area Dimensions	6" W x 2.25" F-B x 6.5" D (20.3 cm x 5.7 cm x 16.5 cm)	
External Dimensions	14.0" W x 13.5" F-B x 14.0" H (35.6 cm x 34.3 cm x 35.6 cm)	
Heater	750 Watts	
Safeties	Over-temperature heater safety & low liquid level alarm	
Reservoir Capacity	1.4 gallons (5.3 liters)	
Bath Reservoir/Lid Material	Stainless steel	
Cabinet Material	Stainless steel and powder coated metal	
Electrical	115V, 60 Hz, 10 Amps	230V, 50/60 Hz, 5 Amps
Shipping Weight	53 lbs. (24.0 kg)	

Specifications are subject to change without notice. 230V, 50/60 Hz unit is CE Marked.