

SPECIFICATION SHEET

VAPORIZED H2O2 DECONTAMINATION CHAMBER

Models: 7800-25-1 / 7800-25-2 / 7800-25-3

SPECIFICATIONS	
Cycle Temperature	45C
Cycle Length	105 Minutes*
Interior Dimensions	32"W x 27"D x 53"H (81 cm x 69 cm x 134 cm)
Exterior Dimensions	36"W x 33"D** x 77"H (90cm x 85cm x196cm)
Interior Construction	2B Finish Solid Stainless Steel (Type 304)
Exterior Construction	Powder Coated Cold Rolled Steel
Work Space	25ft ³ (708L)
Shelf Count	6 Standard
Shelf Construction	Type 304 Stainless Steel Wire, Electropolished
Shelf Dimensions	29.25"W x 26.5"D (74.3cm x 67.3cm)
Shipping Weight	575 lbs(261 kg)/ 875 lbs (397 kg)***

Not for reprocessing of FDA approved medical devices. **Add 2.75" for handle. *Includes export shipping crate.*

ELECTRICAL			
Model	7800-25-1	7800-25-2	7800-25-3
Voltage	115V	208-230V	230V
Frequency	60Hz	60Hz	50Hz
Current	12A	8A	8A

ACCESSORY SHEET

VAPORIZED H₂O₂ DECONTAMINATION CHAMBER

Models: 7800-25-1 / 7800-25-2 / 7800-25-3

ACCESSORIES

Part No.	Description
STER309	Set of 10 Single-use 60ml H ₂ O ₂ Canisters
STER310	Set of 50 Single-use 60ml H ₂ O ₂ Canisters
STER311	Set of 50 Single-use Chemical Indicator (CI) Strips
SHLF346	Additional Stainless Steel Wire Shelf / 2 Shelf Tracks