

CO₂ INCUBATOR / WALLY & WALLYGT

Models: 7420-5 / 7410-5 / 7411-5

Model	7420-5	7410-5	7411-5
Thermal Management	Direct Heat	Gel Jacket	
Temperature Control	±-0.1C at 37°C		
Temperature Uniformity	±0.3°C		
Temperature Range	5°C Above Ambient to 45°C		
Humidity Range	Elevated Up to 90% at 37°C		
Controlled Humidity Range	N/A	Optional Ambient to 90% at 37°C	
CO ₂ Range	Ambient to 20%		
CO ₂ Control	±0.1%		
CO ₂ Sensor	IR		
O ₂ Range	N/A		1-21% O ₂
O ₂ Control	N/A		±-0.1°C
O ₂ Sensor	N/A		Fuel Cell
Contamination Control	Standard ISO5 Chamber HEPA Air Filtration		
H ₂ O ₂ Sterilization	Optional Prep & Module	Standard Prep / Optional Module	
Interior Dimensions	31.3" W x 9.5" D x 26" H (79 cm x 24 cm x 66 cm)		
Exterior Dimensions	36" W x 13.6" D x 36" H (91.5 cm x 34.6 cm x 91.5 cm)		
Interior Construction	Type 304 / Stainless Steel / Electropolished		
Exterior Construction	Powder-Coated Cold Rolled Steel		
Work Space	5 cu.ft. (142 Liters)		
Shelf Count	3 Standard / 13 Maximum		
Shelf Construction	Type 304 / Perforated Stainless Steel / Electropolished		
Shelf Dimensions	31" W x 8.8" D (79 cm x 22 cm)		
Configuration Options	Wall Mount / Benchtop / Dual or Quad Rack		
Unit Installed Weight	200 Lbs. (91 kg.)		
Unit Ship Weight	60 Hz / 227 Lbs. (103 kg.) 50 Hz / 275 Lbs. (125 kg.)*		

*Includes Export Shipping Crate.

CO₂ INCUBATOR / WALLY & WALLYGT

Models: 7420-5 / 7410-5 / 7411-5

Electrical	7420-5-1 7410-5-1 7411-5-1	7420-5-2 7410-5-2 7411-5-2	7420-5-3 7410-5-3 7411-5-3
Volts	115V	208/230V	230V
Hertz	60Hz	60Hz	50Hz
Amps	5A	2A	2A