

CO₂ INCUBATOR / WALLY

Models: 7410-5 / 7411-5

Model	7410-5	7411-5
Temperature Range	5°C Above Ambient to 45°C	
Temperature Control	±0.1°C at 37°C	
Temperature Uniformity	±0.3°C at 37°C	
Humidity Range	Elevated Up to 90% at 37°C	
CO ₂ Range	Ambient to 20%	
CO ₂ Control	±0.1%	
CO ₂ Sensor	IR	
O ₂ Range	N/A	1 to 21% O ₂
O ₂ Control	N/A	±0.1% O ₂
O ₂ Sensor	N/A	Fuel Cell
Interior Dimensions	31.3" W x 9.5" FB x 26" H	
Exterior Dimensions	36" W x 13.6" FB* x 36" H	
Interior Construction	Electropolished Stainless Steel	
Exterior Construction	Powder Coated Cold Rolled Steel	
Work Space	5 cu.ft. (142 liters)	
Shelf Count	3 Standard / 13 Maximum	
Shelf Construction	Perforated Stainless Steel	
Shelf Dimensions	31.3" W x 8.8" FB (79cm x 22cm)	
Configuration Options	Wall Mount / Benchtop Dual, or Quad Rack	
Unit Installed Weight	200 lbs.	
Unit Ship Weight	227 lbs. (60Hz units) 275 lbs. (50Hz units)	

* Add 1.7" FB for handle.

Electrical	
7410-5-1 & 7411-5-1	115V / 60Hz / 5A
7410-5-2 & 7411-5-2	230V / 60Hz / 2A
7410-5-3 & 7411-5-3	230V / 50Hz / 2A