

FINGERPRINT DEVELOPMENT CHAMBERS

Model: 6105

| Model | 6105-2 | 6105-3 |
|-----------------------|---|------------------|
| Temperature Range | 40°C to 100°C | |
| Temperature Control | ±0.5°C | |
| Humidity Range | 40 to 80% RH (up to 80°C) | |
| Humidity Control | ±3.0% RH | |
| Interior Dimensions | 24" W x 16" D x 22" H (61 cm x 41 cm x 56 cm) | |
| Interior Construction | White Epoxy Coated Seam Welded Aluminum | |
| Exterior Dimensions | 29" W x 25" D x 37" H (74 cm x 64 cm x 94 cm) | |
| Exterior Construction | Powder Coated Cold Rolled Steel | |
| Work Space | 5 cu. ft. (140 Liters) | |
| Shelf Work Area | 4.5 sq. ft. (shelf plus bottom) | |
| Electrical | 208/230V, 60 Hz, 15A | 230V, 50 Hz, 15A |
| Shipping Weight | 350 lbs. | 450 lbs.* |

NOTE: Specifications are subject to change without notice.

*Includes export shipping crate.